

BlueCap Washing Buffer based on TRIS without detergents

For use in immunohistochemistry

Product description

| | |
|---------------|--|
| Composition: | TRIS-based buffer with 350 mM NaCl 10x stock solution contains neither Tween nor other detergents |
| pH-value: | pH 7.2 ± 0.2 |
| Preservative: | contains 0.1% ProClin® 300 |
| Storage: | 10x stock solution: shelf life at -20 °C: 1 year - tolerates repeated freezing and thawing cycles shelf life at 2 – 8 °C: 6 months |

Use working solution immediately!

For research use only, not for diagnostic use

| | | |
|----------------|--------|--------------------|
| Available | | |
| Package sizes: | 500 ml | order no. S211.500 |

Instructions for use

This *Washing Buffer* is designed especially for immunohistochemistry. It should be used, whenever detergents can make problems in detection.

Crystals of salt can precipitate after storage at 2 – 8 °C or after freezing. Therefore *Washing Buffer* must be warmed up to room temperature and should be mixed thoroughly before preparing the working solution. This leads to dissolving of salt after shaking. Stock solution is diluted 1:10 with salt free water to get the working solution.

Use working dilution immediately.

ProClin is a registered trademark of Rohm and Haas Company.