

Try Before You Buy · EasySort™ Cell Isolation Kits

Free 10 assays

Complete Cell Isolation Solution


Up to \$852 in Benefits

Voucher rewards

Promotion: May 1st-Aug 31st, 2026

Why Subset Isolation Matters

- ★ Functional differences drive outcomes
- ★ Purity affects data quality
- ★ Exhaustion/memory studies require precision



Our Kits Deliver

- High purity & recovery
- Fast workflow (~14-18 min)
- Gentle on cells
- Ready for functional assays

Isolate the exact T cell subset you need — naïve, effector, memory, or regulatory

FREE Claim Your Free EasySort™ Kit

Step 1



Free Trial

10 assays + magnetic stand EC001 (\$852 value)

Submit Feedback

\$50 voucher per report

Step 2



Step 3



Upgrade

Purchase ≥500 assays → keep the stand **FREE**

Apply for trial pack through your local distributor

Choose From!



Product Name	Cat. No.	Size
EasySort™ Human Naïve Pan T Cell Isolation Kit	MIH006N	100/200 Assays
EasySort™ Human Naïve CD4+ T Cell Isolation Kit	MIH007N	100/200 Assays
EasySort™ Human Naïve CD8+ T Cell Isolation Kit	MIH008N	100/200 Assays
EasySort™ Human CD3+ T Cell Isolation Kit	MIH001N	100/200 Assays
EasySort™ Human CD4+ T Cell Isolation Kit	MIH002N	100/200 Assays
EasySort™ Human CD8+ T Cell Isolation Kit	MIH003N	100/200 Assays
EasySort™ Human B Cell Isolation Kit	MIH004N	100/200 Assays
EasySort™ Human Memory CD4+ T Cell Isolation Kit	MIH009N	100/200 Assays
EasySort™ Human Memory CD8+ T Cell Isolation Kit	MIH010N	100/200 Assays
EasySort™ Mouse Pan-Naïve T Cell Isolation Kit	MIM006N	100/200 Assays
EasySort™ Mouse Naïve CD4+ T Cell Isolation Kit	MIM007N	100/200 Assays
EasySort™ Mouse Naïve CD8+ T Cell Isolation Kit	MIM008N	100/200 Assays
EasySort™ Mouse CD3+ T Cell Isolation Kit	MIM001N	100/200 Assays
EasySort™ Mouse CD4+ T Cell Isolation Kit	MIM002N	100/200 Assays
EasySort™ Mouse CD8+ T Cell Isolation Kit	MIM003N	100/200 Assays
EasySort™ Mouse B Cell Isolation Kit	MIM004N	100/200 Assays
EasySort™ Mouse NK Cell Isolation Kit	MIM005N	100/200 Assays

Validated Performance

Human CD3+ T Cells
60.9% → 96.7% purity

Mouse CD4+ T Cells
25.5% → 95.9% purity

