

買Kits就送 星巴克禮品卷



ELISA Kits 優勢



- ★ 快速且容易使用：Quickey Pro™ ELISA 只需1.5小時
- ★ 低樣品量需求：Mini Sample ELISA只需25ul
- ★ 超級靈敏



ELISA 試劑盒產品比較表

項目	傳統 ELISA 試劑盒	Quickey Pro ELISA	Mini Sample ELISA 試劑盒	高靈敏ELISA試劑盒
主要優勢	一般常見的蛋白	最少時間	極少量樣本檢測	高靈敏度
偵測方式	Sandwich/Competitive	Sandwich/Competitive	Sandwich	Sandwich
樣品體積	50μL /100 μL	50 μL	25 μL	100 μL
操做時間	2.5h/3.5h	1.5h	3.5h	3.5h
所需步驟	2步/3步	1步	3步	3步
檢測靈敏度	pg/mL	pg/mL	pg/mL	fg/mL
儲存方式	2-8度/12個月	2-8度/6個月	2-8度/12個月	2-8度/12個月



I Quickey Pro™ ELISA試劑盒特色

產品特性

- ⚡ 操作時間短：完成實驗僅需90分鐘
- 🕒 操作簡單：簡化流程僅需1步驟清洗/添加樣品
- ↓ 僅需少量樣本：每次測試僅需50μL

應用場景

- 1 適用於你的實驗時間緊迫的時候
- 2 適合樣本有限的時候，可以檢測多個指標
- 3 能檢測特殊物種，包括猴子、豬、兔子、雞、牛等等

I Mini Sample ELISA Kits試劑盒特色

產品特性

- ↓ 所需樣品量最少：每次測試僅需25ul
- 🔄 優化過的生產流程：提供更好的靈敏度跟提升效能
- 🔬 研究驅動創新：持續開發新產品符合研究趨勢

應用場景

- 1 適合週期或階段性研究
- 2 可以在一隻動物內偵測各種多樣化指標
- 3 適用於稀有動物或高價值樣本

I High Sensitivity ELISA Kits

產品特性

- ↑ 超高靈敏度：最低可檢測至0.2pg/mL
- 📈 提高檢測率：在normal serum/plasma中可達30-100%
- 🧪 多樣本類型驗證：適用於一般分泌形樣本和細模式樣本

應用場景

- 1 適合低濃度樣本需要超高靈敏度時
- 2 專為偵測細胞因子設計，例如IL、IFN、TNF等等

I High-Impact Publications Featuring Our Products

- ① A Novel Antidiuretic Hormone Governs Tumour-Induced Renal Dysfunction. (*Nature*, 2023)

Products:

- cGMP ELISA Kit (E-EL-0083)
- Mouse ADH ELISA Kit (E-EL-M0106)
- Mouse ALB ELISA Kit (E-EL-M3032)
- Mouse BMG/β2-MG ELISA Kit (E-EL-M2411)

- ② A Gut-Derived Hormone Regulates Cholesterol Metabolism. (*Cell*, 2024)

Products:

- Mouse ApoB ELISA Kit (E-EL-M3017)



TEL : 02-8866-3366
FAX : 02-2833-1789
Mail : info@abscience.com.tw



Line ID:@asiabio